

IMPORTANT



Additional Maintenance for Forensic Buddy®

- The inside of the red housing should be wiped out periodically to prevent the buildup of green powder. We recommend cleaning the unit after 20 rounds. This frequency should be increased if you are firing large volumes of traditionally ‘dirty’ calibers.
- After 50 – 60 rounds, we recommend washing the media in a bucket of D-Lead solution. This will remove any grit / powder that has been deposited on the media. After washing the media in the D-Lead, it should be rinsed in clean water and allowed to dry. The D-Lead solution can be reused many times.

3/3/2010